

Accessories, apps, & hearing aid compatibility.

Hearing aid	Charging stations		Performance levels			CROS/BICROS				App	Remote control	Programming	
	Portable	Li Charger	Premium	Advanced	Standard	CROS RIC LI 7	CROS RIC 6	CROS RIC LI 6	CROS RIC 5			TruHearing app	miniPocket
New TH RIC LI 7	•	•	•	•	-	**	-	-	-	•	-	•	-
New TH RIC LI 7 (T)	•	•	•	•	-	**	-	-	-	•	-	•	-
New TH IF 7	•	-	•	-	-	**	-	-	-	•	•	-	•
TH RIC LI 6	-	•	•	•	-	-	-	**	-	•	•	•	-
TH RIC LI 6 (T)	-	•	•	•	-	-	**	**	*	•	•	•	-
TH RIC 6	-	-	•	•	•	-	**	-	*	•	•	•	-
TH SR LI 6	•	•	•	-	-	-	**	-	-	•	•	•	-
TH BTE LI 5	-	•	•	•	-	-	-	-	•	•	•	•	-
TH Power BTE LI 5	-	•	•	•	-	-	-	-	•	•	•	•	-
SP BTE LI 5	-	•	•	•	-	-	-	-	•	•	•	•	-
TH 6 LI ITE/HS/ITC	-	•	•	-	-	-	**	**	-	•	•	•	-
TH IIC 19	-	-	-	-	•	-	-	-	-	•	•	-	•
TH CIC 19	-	-	-	-	•	-	-	-	-	•	•	-	•
TH 19 ITE/HS/ITC	-	-	-	-	•	-	-	-	-	•	•	•	•
TH 19 STD BTE	-	-	-	-	•	-	-	-	-	•	•	•	-
TH 19 Power BTE	-	-	-	-	•	-	-	-	-	•	•	•	-

* Recommended compatibility for ease of patient use † Only compatible with Bluetooth configurations ** CROS compatible version available



(RIC) Receiver-in-canal



(BTE) Behind the ear



(CIC/IIC) Custom



Questions on products?

Contact us at the dedicated TruHearing support line
1-800-350-6147 or email TruHearing@signiausa.com

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Signia is under license. Other trademarks and trade names are those of their respective owners.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc. iOS is a trademark or registered trademark of Apple Inc. in the U.S. and other countries and is used under license.

TruHearing® is a registered trademark of TruHearing, Inc.

Copyright © 2024 Signia is a registered trademark of WSAUD A/S. SIV-10240588



TruHearing®

Portfolio overview Spring 2024



Finally, a hearing aid designed for real conversations.

Meet the world's first hearing aid with a multi-stream architecture capable of accurately pinpointing multiple moving speakers in real-time conversation.

We developed the TruHearing 7 platform by leveraging our deep understanding of dynamic conversations, including head movements of the wearer and changing positions of speakers. This groundbreaking technology ensures uninterrupted conversation engagement in noise, addressing a critical challenge many hearing aid wearers face.

Future-proof connectivity

Get future-proofed connectivity with TruHearing 7. Our devices are ready for the upcoming Bluetooth® LE Audio standard. Give wearers the chance to stay connected with the future of hearing care.

- ASHA (Audio Streaming for Hearing Aids)*
- Made for iPhone®
- LE Audio ready* via future firmware update
- Bluetooth® connectivity via StreamLine™ Mic*



The rechargeable and discreet TruHearing 7 Instant Fit provides power to hear all day.

The TruHearing 7 Instant Fit style gives patients the convenience of a rechargeable battery for all day listening, plus the accompanying case charges wirelessly, adding an extra level of hassle-free convenience. Take advantage of up to 28 hours of use on a single charge, with a quick 30-minute quick-charge option. The convenient charging case gives up to four charges on the go.

CROS TruHearing 7 Instant Fit is also now available, giving you even more flexibility in providing your patients with personalized hearing solutions.

*TruHearing 7 Instant Fit limited to use of TruHearing app only.

Ready-to-wear & customs



TruHearing 7 Instant Fit Rechargeable




TruHearing CIC/IIC 19




TruHearing 6 LI ITE/HS/ITC Rechargeable



TruHearing CA/HS/ITE 19



TruHearing ITE 19



Optional Telecoil

Bluetooth

New rechargeable RICs



TruHearing RIC LI 7



TruHearing RIC LI 7 (T)

Rechargeable RICs & BTEs



TruHearing RIC LI 6 (T)



TruHearing RIC LI 6



TruHearing SR LI 6 (SLIM RIC)




TruHearing BTE LI 5



TruHearing Power BTE LI 5



TruHearing SP BTE LI 5



For your patients with tinnitus, three strategies are available to provide relief.



World's first Tinnitus Notch Therapy



Static Noise



Ocean Waves

The science of keeping up with conversations.

The revolutionary TruHearing 7 platform with Dynamic Conversation Tracking* analyzes, augments, and adapts to the dynamic flow of conversations. This lets wearers enjoy a seamless and enriched experience.

Analyze

The platform processes 192,000 data points every second to identify speakers and analyze conversation dynamics, authentically reflecting real-life conversations.

Augment

Uses data gathered to control multiple focus streams that lock on to and dynamically highlight each speaker's voice, even amidst noise. Meanwhile, (Re)Mix technology processes background sound to ensure a comfortable, immersive conversational experience without distractions.

Adapt

Our multi-stream architecture updates 1,000 times per second, and focus streams merge into a dynamic auditory space that adapts to evolving conversations. This allows for real-time tracking and enhancement of multiple moving speakers, ensuring wearers can fully engage in group conversations.

*Dynamic Conversation Tracking not available in all styles.

