

TH 7 Products

Features & Accessories

Data Sheet

TH 7 RIC LI

TH Premium 7 RIC LI TH Advanced 7 RIC LI TH Demo 7 RIC LI

TH 7 RIC LIT

TH Premium 7 RIC LI T TH Advanced 7 RIC LI T

TH 7 RIC LI T BC

TH Premium 7 RIC LI T BC TH Advanced 7 RIC LI T BC TH Standard 7 RIC LI T BC

TH 7 IF LI

TH Premium 7 IF LI

TH 7 SR LI

TH Premium 7 SR LI

TH 7 Canal/ITE/CIC/IIC

TH Premium 7 Canal TH Advanced 7 Canal TH Standard 7 Canal

TH Premium 7 ITE TH Advanced 7 ITE TH Standard 7 ITE

TH Premium 7 CIC TH Advanced 7 CIC TH Standard 7 CIC

TH Premium 7 IIC TH Advanced 7 IIC TH Standard 7 IIC

TH 7 RIC LI TH 7 RIC LI T TH 7 RIC LI T BC

Features

	Premium	Advanced	Standard
Conversation Management			
RealTime Conversation enhancement 1)			_
Conversation dynamics analyzer	•	_	_
Conversation source analyzer	•	•	_
Dynamic Focus streams	3	1	_
Conversation booster			_
Augmented Focus	•	•	•
Focus streams	1	1	1
Surrounding streams	1	1	1
Dynamic Soundscape Processing 3.0			
Spatial SpeechFocus 1)	•	_	_
SpeechFocus	•	•	_
Own Voice Processing 2.0 (OVP 2.0) 1)	•	•	_
Motion sensor	•	•	_
Binaural OneMic Directionality 2.0 1)	_	_	_
Binaural directionality 1)	•	•	_
Adaptive directionality	•	•	•
Automatic directionality	•	•	•
TruEar 360	•	•	•
SpeechStabilizer	•	•	•
Speech and noise management	•	•	•
e2e wireless 4.0 ¹⁾	•	•	•
Listening Comfort		'	
Extended bandwidth	•	_	_
Auto EchoShield	•	_	_
Reverberant room program	•	•	_
SoundSmoothing			
eWindScreen			
Extended dynamic range	•	•	•
Feedback cancellation	•	•	•
Wearer Empowerment		'	
Spatial Configurator	•	•	_
My WellBeing	•	•	<u>—</u>
TruHearing App (iOS, Android)	•	•	•
Adaptive streaming volume	•	•	•
TeleCare	•	•	•

Note: TH Demo 7 RIC LI can be set to Premium only.

TH 7 RIC LI TH 7 RIC LI T TH 7 RIC LI T BC

Features

	Premium	Advanced	Standard
Connectivity			
Direct streaming for iOS (MFi) & Android (ASHA)	•	•	•
HandsFree for iOS	•	•	•
CallControl for iOS	•	•	•
Telecoil 2)	•	•	•
StreamLine Mic compatibility	•	•	_
StreamLine TV compatibility	•	•	
miniPocket compatibility	•	•	•
Tinnitus			
Ocean Waves Therapy signal (presets)	4	4	4
Static Therapy signal (channels / presets)	20 / 5	18 / 5	12 / 5
Notch Therapy	•	•	•
Tinnitus noise therapy manager	•	•	•
Fitting Optimization			
Processing channels / Gain & MPO handles	48 / 20	36 / 18	24 / 12
HD Music (presets)	3	3	1
Hearing programs	6	6	4
CROS/BiCROS compatibility	•	•	_
Frequency compression	•	•	•
Smart Optimizer and Data Logging	•	•	•
Acclimatization manager	•	•	•
InSituGram	•	•	•
AutoFit	•	•	•
 available for bilateral, wireless fittings TH 7 RIC LI T & TH 7 RIC LI T BC only 	availablenot available	highest featur	re performance

Note: TH Demo 7 RIC LI can be set to Premium only.

TH 7 ITE/Canal⁵⁾

	Premium	Advanced	Standard
Conversation Management			
RealTime Conversation Enhancement 1)			_
Conversation source analyzer	•	•	_
Conversation dynamics analyzer	•	_	_
Dynamic Focus streams	3	1	_
Conversation booster			_
Augmented Focus	•	•	•
Focus streams	1	1	1
Surrounding streams	1	1	1
Dynamic Soundscape Processing 3.0			
Spatial SpeechFocus 1)	•	_	_
SpeechFocus	•	•	_
Own Voice Processing 2.0 (OVP 2.0) 1)		_	_
Motion sensor		_	_
Binaural OneMic Directionality 2.0 1)	_	_	_
Binaural directionality 1)	•	•	_
Adaptive directionality	•	•	•
Automatic directionality	•	•	•
TruEar 360	•	•	•
SpeechStabilizer	•	•	•
Speech and noise management	•	•	•
e2e wireless 4.0 ¹⁾⁴⁾	•	•	•
Listening Comfort		'	'
Feedback cancellation	•	•	•
eWindScreen			
SoundSmoothing			
Auto EchoShield	•	_	_
Reverberant room program	•	•	
Extended dynamic range	•	•	•
Extended bandwidth	•	_	_
Wearer Empowerment			
Spatial Configurator	•	•	_
My WellBeing		_	_
TruHearing App (iOS, Android)	•	•	•
Adaptive streaming volume		_	_
TeleCare 3)	•	•	•
		·	

¹⁾ available for bilateral fittings with e2e wireless

²⁾ for remote control adjustments

³⁾ for remote services TeleCare Dashboard and CareChat

⁴⁾ for devices with e2e wireless link. Not available for non-wireless TH 7 ITE/Canal/CIC/IIC

⁵⁾ two-microphone versions

available

⁻ not available

optional

TH 7 ITE/Canal⁵⁾

	Premium	Advanced	Standard
Connectivity			
Direct streaming for iOS (MFi) & Android (ASHA)	_	_	_
HandsFree for iOS	_	_	_
CallControl for iOS	_	_	_
Telecoil	0	0	0
StreamLine Mic compatibility	_	_	_
StreamLine TV compatibility	_	_	_
miniPocket compatibility	•	•	•
Tinnitus		'	'
Ocean Waves Therapy signal (presets)	4	4	4
Static Therapy signal (channels / presets)	20 / 5	18 / 5	12 / 1
Notch Therapy	•	•	•
Tinnitus noise therapy manager	•	•	•
Fitting Optimization		'	
Processing channels / Gain & MPO handles	48 / 20	36 / 18	24 / 12
HD Music (presets)	3	3	1
Hearing programs	6	6	4
CROS/BiCROS compatibility 4)	•	•	•
Frequency compression	•	•	•
Smart Optimizer and Data Logging	•	•	•
Acclimatization manager	•	•	•
InSituGram	•	•	•
AutoFit	•	•	•
		·	

¹⁾ available for bilateral fittings with e2e wireless link

²⁾ for remote control adjustments

³⁾ for remote services TeleCare Dashboard and

⁴⁾ for devices with e2e wireless link. Not available for non-wireless TH 7 ITE/Canal/CIC/IIC

⁵⁾ two-microphone versions

available

not available

optional

TH 7 IF LI TH 7 Canal⁵⁾/CIC/IIC

	Premium	Advanced 5)	Standard 5)
Conversation Management			
RealTime Conversation enhancement 1)	_	_	_
Conversation dynamics analyzer	<u> </u>	_	_
Conversation source analyzer		_	_
Dynamic Focus streams		_	_
Conversation booster		_	_
Augmented Focus	_	_	_
Focus streams	_	_	_
Surrounding streams	_	_	_
Dynamic Soundscape Processing 3.0			
Spatial SpeechFocus ¹⁾		_	_
SpeechFocus	_	_	_
Own Voice Processing 2.0 (OVP 2.0) ¹⁾	_	_	_
Motion sensor	_	_	_
Binaural OneMic Directionality 2.01)	•	•	_
Binaural directionality 1)	_	_	_
Adaptive directionality	_	_	_
Automatic directionality	<u> </u>	_	_
TruEar 360	<u> </u>	_	_
SpeechStabilizer	•	•	•
Speech and noise management	•	•	•
e2e wireless 4.0 ¹⁾⁴⁾	•	•	•
Listening Comfort		'	
Extended bandwidth	•	•	•
Auto EchoShield			
Reverberant room program			
SoundSmoothing	•	_	_
eWindScreen	•	•	_
Extended dynamic range	•	•	•
Feedback cancellation	•	_	_
Wearer Empowerment		'	
Spatial Configurator	-	_	_
My WellBeing	<u> </u>	_	_
TruHearing App (iOS, Android) ²⁾	•	•	•
Adaptive streaming volume		_	_
TeleCare 3)	•	•	•
1) available for bilateral, wireless fittings	available	highest featur	re nerformance

¹⁾ available for bilateral, wireless fittings

²⁾ for remote control and support via TeleCare

³⁾ for remote services TeleCare Dashboard and CareChat

⁴⁾ for devices with e2e wireless link. Not available for TH 7 CIC/IIC non-wireless

⁵⁾ one-microphone versions

available

not available

optional

highest feature performance

TH 7 IF LI TH 7 Canal⁵⁾/CIC/IIC

	Premium	Advanced 5)	Standard 5)
Connectivity			
Direct streaming for iOS (MFi) & Android (ASHA)	_	_	_
HandsFree for iOS	_	_	
CallControl for iOS	_	_	_
Telecoil	0	0	0
StreamLine Mic compatibility	_	_	_
StreamLine TV compatibility		_	-
miniPocket compatibility	•	•	•
Tinnitus			
Ocean Waves Therapy signal (presets)	4	4	4
Static Therapy signal (channels / presets)	20 / 5	18 / 5	12 / 5
Notch Therapy	•	•	•
Tinnitus noise therapy manager	•	•	•
Fitting Optimization			
Processing channels / Gain & MPO handles	48 / 20	36 / 18	24 / 12
HD Music (presets)	3	3	1
Hearing programs	6	6	4
CROS/BiCROS compatibility ⁴⁾	•	•	•
Frequency compression	•	•	•
Smart Optimizer and Data Logging	•	•	•
Acclimatization manager	•	•	•
InSituGram	•	•	•
AutoFit	•	•	•
1) available for bilateral, wireless fittings	available	highest featur	e performance

¹⁾ available for bilateral, wireless fittings

²⁾ for remote control and support via TeleCare

³⁾ for remote services TeleCare Dashboard and CareChat

⁴⁾ for devices with e2e wireless link. Not available for TH 7 CIC/IIC non-wireless

⁵⁾ one-microphone versions

available

not available

optional

TH 7 SR LI

	Premium
Conversation Management	
RealTime Conversation enhancement 1)	
Conversation dynamics analyzer	•
Conversation source analyzer	•
Dynamic Focus streams	3
Conversation booster	
Augmented Focus	•
Focus streams	1
Surrounding streams	1
Dynamic Soundscape Processing 3.0	
Spatial SpeechFocus 1)	•
SpeechFocus	•
Own Voice Processing 2.0 (OVP 2.0) 1)	•
Motion sensor	•
Binaural OneMic Directionality 2.0 1)	_
Binaural directionality 1)	•
Adaptive directionality	•
Automatic directionality	•
TruEar 360	•
SpeechStabilizer	•
Speech and noise management	•
e2e wireless 4.0 ¹⁾	•
Listening Comfort	
Extended bandwidth	•
Auto EchoShield	•
Reverberant room program	•
SoundSmoothing	
eWindScreen	
Extended dynamic range	•
Feedback cancellation	•
Wearer Empowerment	
Spatial Configurator	•
My WellBeing	•
TruHearing App (iOS, Android) ²⁾	•
Adaptive streaming volume	•
TeleCare 3)	•
1) available for bilateral, wireless fittings 2) for remote control and support via TeleCare	available not available

²⁾ for remote control and support via TeleCare

³⁾ for remote services TeleCare Dashboard and CareChat

⁻ not available

TH 7 SR LI

| Features

	Premium
Connectivity	
Direct streaming for iOS (MFi) & Android (ASHA)	•
HandsFree for iOS	•
CallControl for iOS	•
Telecoil	_
StreamLine Mic compatibility	•
StreamLine TV compatibility	•
miniPocket compatibility	•
Tinnitus	
Ocean Waves Therapy signal (presets)	4
Static Therapy signal (channels / presets)	20 / 5
Notch Therapy	•
Tinnitus noise therapy manager	•
Fitting Optimization	
Processing channels / Gain & MPO handles	48 / 20
HD Music (presets)	3
Hearing programs	6
CROS/BiCROS compatibility	•
Frequency compression	•
Smart Optimizer and Data Logging	•
Acclimatization manager	•
InSituGram	•
AutoFit	•
1) available for bilateral, wireless fittings	available

²⁾ for remote control and support via TeleCare

³⁾ for remote services TeleCare Dashboard and CareChat

⁻ not available

Configuration and Accessories

		TH 7	TH 7	TH 7
Performance level	Premium Advanced, Standard	RIC LI	RIC LI T	RIC LI T BC
Instrument and compo	nents			
Ingress Protection rating		IP68	IP68	IP68
Battery type/size		Lithium-ion	Lithium-ion	Lithium-ion
Programming accessor	ries			
Noahlink Wireless / Noah	nlink Wireless 2	_	A	_
ConnexxAir		_	_	_
Programming adapter / c	cable	_	_	_
Accessories				
miniPocket		0	0	0
StreamLine TV		0	0	_
TV Sound		<u> </u>	_	0
StreamLine Mic		0	0	_
TruHearing App (iOS, Ar	ndroid)	0	0	0
Charger			A	A
Multi Charger		0	0	0
Portable Charger R	IC	0	0	0
Charging Station R		0	0	0
Charging+ Station F	?	0	0	_
IF CIC Travel Charg	ger	_	_	_
CROS Pure C&G IX		0	0	_
CROS Silk C&G IX		_	_	_
CROS Pure C&G BCT IX	<	_	_	0
		mandatory	 not available 	

optional

Configuration and Accessories

		TH 7
Performance level	Premium Advanced	IF LI
Instrument and components		
Ingress Protection rating		IP68
Battery type/size		Lithium-ion
Programming accessories		
Noahlink Wireless / Noahlink Wirel	ess 2	_
ConnexxAir		A
Programming adapter / cable		_
Accessories		
miniPocket		0
StreamLine TV		_
StreamLine Mic		_
TruHearing App (iOS, Android)		0
Charger		A
Multi Charger		_
Portable Charger RIC		_
Charging Station R		_
Charging+ Station R		_
IF CIC Travel Charger		A
CROS Pure C&G IX		_
CROS Silk C&G IX		0
		▲ mandatory

optional

Configuration and Accessories

		TH 7	TH 7	TH 7
Performance level	Premium Advanced, Standard	Canal/ITE	CIC/IIC	SR LI
Instrument and compor	nents			
Ingress Protection rating		_	_	IP68
Battery type/size		312	10	Lithium-ion
Programming accessor	ies			
Noahlink Wireless / Noah	link Wireless 2	_	_	
ConnexxAir		○ 2)	<u> </u>	_
Programming adapter / ca	able	0	Flex connector (T5 FP) / Adapter (I3 FP)	_
Accessories				
miniPocket		0	0	0
StreamLine TV		_	_	0
StreamLine Mic		_	_	0
TruHearing App (iOS, An	droid)	0	0	0
Charger		_	_	A
Universal Cable Cha	arger	_	_	_
Travel Charger RIC		_	_	_
Charging Station R		_	_	_
Charging+ Station R		_	_	_
SR Travel Charger			_	A
CROS Pure C&G IX		_		_
CROS Silk C&G IX		<u> </u>	<u> </u>	_
CROS Styletto IX		_	_	0
1) TH 7 CIC/IIC wireless 2) e2e wireless link only	only	mandatoryoptional		

not available